# Easy Line PartnerMic

This easy-to-use lapel worn microphone is ideal for one-to-one conversations in noisy environments. It can be attached anywhere close to the speakers voice regardless of directionality or placement. With integrated Sonova AirStream™ technology, clients can benefit from clear speech streamed to compatible hearing aids.

# Technical data

Dimensions (LxWxH)	55 x 30 x 18 mm (2.14 x 1.2 x 0.7 ")			
Weight	18.2 g (0.04 lbs)			
Color	Grey			
Operating conditions	0° to +45° Celsius (+32° to +113° Fahrenheit). Relative humidity of <90% (non-condensing)			
Transport and storage conditions	-20° to +60° Celsius (-4° to +140° Fahrenheit). Relative humidity of <95% for a long period of time			




### Features

- Out of the box pairing
- Always directional regardless of position
- Sonova AirStream technology

### Device description

1	Microphone openings
2	Charging input
3	Clip
4	Pairing button
5	Power switch On/Off
6	Indicator light (LED)







Battery	ir	nfo	rm	atio	n

Dattery IIIIOIIIIation					
Туре	Lithium Polymer 3.7 VDC. Embedded, not removable				
Capacity	230 mAh (typ.)				
Operating time	Up to 6 hours streaming				
Charging time	2 hours (typ.)				
Power supply					
Туре	DC, stabilized				
Secondary voltage	5 VDC				
Max. current	1 A				
Plug	USB-C				
Important notice	Only use the original power supply				
Audio information					
Audio quality	120 Hz - 7.3 kHz				
Radio information					
Streaming technology	AirStream technology				
Frequency range	2.4 GHz ISM band				
RF power	20 mW max. (conducted)				
Transmission range	Up to 25 meters / 80 feet				
Max. connected devices	1 pair of hearing aids				
Standards					
Europe	EN 300 328; EN 301 489; IEC/EN 60601-1-2;				
	IEC/EN 62479; EN/IEC 62368-1; IEC 62133				
Canada	RSS-247; RSS-102				
Japan	ARIB-T66				
USA	CFR 47, part 2.1093; KDB 447498				

## Compatibility note

Easy Line PartnerMic is only compatible with Sonova direct streaming hearing aids.





And importer for European Union: Sonova Deutschland GmbH Max-Eyth-Str. 20 70736 Fellbach-Oeffingen, Germany

